

Product Datasheet

Human Extracellular Superoxide Dismutase (Cu-Zn) ELISA Kit (orb1213535)

Description	Human Extracellular Superoxide Dismutase (Cu-Zn), SOD3 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	SOD3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,851
Uniprot ID	P08294
Solubility (25°C)	960ng/L
Sensitivity	1.45ng/L
Entrez	6649
Expiration Date	Please enquire.