

Product Datasheet

Human Keratin ELISA Kit (orb1213538)

Description	Human Keratin, Type I Cytoskeletal 16, KRT16 ELISA Kit
Reactivity	Human
Range	0.05-9ng/ml
Tested Applications	ELISA
Target	KRT16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,268
Uniprot ID	P08779
Solubility (25°C)	9.6ng/ml
Sensitivity	0.023ng/ml
Entrez	3868
Expiration Date	Please enquire.