

Product Datasheet

Human Flotillin-1 ELISA Kit (orb1213540)

Description	Human Flotillin-1, FLOT1 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	FLOT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,355
Uniprot ID	075955
Solubility (25°C)	1600ng/L
Sensitivity	3.09ng/L
Entrez	10211
Expiration Date	Please enquire.