

Product Datasheet

Human Netrin Receptor Unc5D ELISA Kit (orb1213554)

Description	Human Netrin Receptor Unc5D, UNC5D ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	UNC5D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	105,880
Uniprot ID	Q6UXZ4
Solubility (25°C)	48ng/ml
Sensitivity	0.046ng/ml
Entrez	137970
Expiration Date	Please enquire.