

Product Datasheet

Human Rho-associated Protein Kinase 2 ELISA Kit (orb1213557)

Description	Human Rho-associated Protein Kinase 2, ROCK2 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	ROCK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	160,900
Uniprot ID	075116
Solubility (25°C)	64ng/ml
Sensitivity	0.08ng/ml
Entrez	9475
Expiration Date	Please enquire.