

Product Datasheet

Human Hyaluronan Synthase 1 ELISA Kit (orb1213560)

Description	Human Hyaluronan Synthase 1, HAS1 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	HAS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	64,832
Uniprot ID	Q92839
Solubility (25°C)	9600ng/L
Sensitivity	15.55ng/L
Entrez	3036
Expiration Date	Please enquire.