

Product Datasheet

Human Bleomycin Hydrolase ELISA Kit (orb1213565)

Description	Human Bleomycin Hydrolase, BLMH ELISA Kit
Reactivity	Human
Range	5-1800ng/L
Tested Applications	ELISA
Target	BLMH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,562
Uniprot ID	Q13867
Solubility (25°C)	1920ng/L
Sensitivity	2.62ng/L
Entrez	642
Expiration Date	Please enquire.