

Product Datasheet

Human Nadphoxidase 4 ELISA Kit (orb1213567)

Description	Human Nadphoxidase 4, NOX4 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	NOX4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	66,932
Uniprot ID	Q9NPH5
Solubility (25°C)	1600ng/L
Sensitivity	2.36ng/L
Entrez	50507
Expiration Date	Please enquire.