

Product Datasheet

Human Luteinizing Hormone Receptor ELISA Kit (orb1213578)

Description	Human Luteinizing Hormone Receptor, LHR ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	TRHR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,085
Uniprot ID	P34981
Solubility (25°C)	32ng/ml
Sensitivity	0.029ng/ml
Entrez	7201
Expiration Date	Please enquire.