

Product Datasheet

Human Toll-like Receptor 6 ELISA Kit (orb1213581)

Description	Human Toll-like Receptor 6, TLR6 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	TLR6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	91,880
Uniprot ID	Q9Y2C9
Solubility (25°C)	12.8ng/ml
Sensitivity	0.01ng/ml
Entrez	10333
Expiration Date	Please enquire.