

Product Datasheet

Human Beta-Actin ELISA Kit (orb1213585)

Description	Human Beta-Actin, B-ACTIN ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	TCTB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,956
Uniprot ID	Q96PL2
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	6975
Expiration Date	Please enquire.