

Product Datasheet

Human Cystic Fibrosis Transmembrane Conductance Regulator ELISA Kit (orb1213586)

Description	Human Cystic Fibrosis Transmembrane Conductance Regulator, CFTR ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CFTR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	168,142
Uniprot ID	P13569
Solubility (25°C)	4000ng/L
Sensitivity	9.28ng/L
Entrez	1080
Expiration Date	Please enquire.