

Product Datasheet

Human Cystatin-A ELISA Kit (orb1213587)

Description	Human Cystatin-A, CSTA ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	CSTA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,006
Uniprot ID	P01040
Solubility (25°C)	32ng/ml
Sensitivity	0.028ng/ml
Entrez	1475
Expiration Date	Please enquire.