

Product Datasheet

Human Ferritin Light Chain ELISA Kit (orb1213589)

Description	Human Ferritin Light Chain, FTL ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	FTL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,020
Uniprot ID	P02792
Solubility (25°C)	96ng/ml
Sensitivity	0.031ng/ml
Entrez	2512
Expiration Date	Please enquire.