

Product Datasheet

Human Protein-arginine Deiminase Type-4 ELISA Kit (orb1213594)

Description	Human Protein-arginine Deiminase Type-4, PADI4 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	PADI4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	74,079
Uniprot ID	Q9UM07
Solubility (25°C)	24ng/ml
Sensitivity	0.021ng/ml
Entrez	23569
Expiration Date	Please enquire.