

Product Datasheet

Human Defensin-6 ELISA Kit (orb1213597)

Description	Human Defensin-6, DEFA6 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	DEFA6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,975
Uniprot ID	Q01524
Solubility (25°C)	24ng/ml
Sensitivity	0.024ng/ml
Entrez	1671
Expiration Date	Please enquire.