



Product Datasheet

APP antibody (orb12136)


Description	Rabbit polyclonal antibody to APP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic peptide of human APP(NP_000475.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100-140KDa
Uniprot ID	P05067
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of rat spleen tissue.



Immunohistochemical staining of rat brain tissue.