

Product Datasheet

Human Interferon Kappa ELISA Kit (orb1213603)

Description	Human Interferon Kappa, IFNK ELISA Kit
Reactivity	Human
Range	1-300ng/L
Tested Applications	ELISA
Target	IFNK
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,218
Uniprot ID	Q9P0W0
Solubility (25°C)	320ng/L
Sensitivity	0.53ng/L
Entrez	56832
Expiration Date	Please enquire.