

Product Datasheet

Human Meteorin-like Protein ELISA Kit (orb1213623)

Description	Human Meteorin-like Protein, METRNL ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	METRNL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,398
Uniprot ID	Q641Q3
Solubility (25°C)	16ng/ml
Sensitivity	0.023ng/ml
Entrez	284207
Expiration Date	Please enquire.