

Product Datasheet

Human Laminin Subunit Alpha-2 ELISA Kit (orb1213624)

Description	Human Laminin Subunit Alpha-2, LAMA2 ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	LAMA2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	343,905
Uniprot ID	P24043
Solubility (25°C)	128ng/ml
Sensitivity	0.29ng/ml
Entrez	3908
Expiration Date	Please enquire.