

Product Datasheet

Human Glutathione Peroxidase 3 ELISA Kit (orb1213627)

Description	Human Glutathione Peroxidase 3, GPX3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	GPX3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,552
Uniprot ID	P22352
Solubility (25°C)	24ng/mL
Sensitivity	0.026ng/mL
Entrez	2878
Expiration Date	Please enquire.