

Product Datasheet

Human Fibroblast Growth Factor 8 ELISA Kit (orb1213630)

Description	Human Fibroblast Growth Factor 8, FGF-8 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	FGF8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,525
Uniprot ID	P55075
Solubility (25°C)	2400ng/L
Sensitivity	4.94ng/L
Entrez	2253
Expiration Date	Please enquire.