

Product Datasheet

Human Secreted Frizzled Related Protein 3 ELISA Kit (orb1213635)

Description	Human Secreted Frizzled Related Protein 3, SFRP3 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	FRZB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,254
Uniprot ID	Q92765
Solubility (25°C)	40ng/ml
Sensitivity	0.021ng/ml
Entrez	2487
Expiration Date	Please enquire.