

Product Datasheet

Human Neurogranin ELISA Kit (orb1213637)

Description	Human Neurogranin, NRGN ELISA Kit
Reactivity	Human
Range	5-1500ng/ml
Tested Applications	ELISA
Target	NRGN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,618
Uniprot ID	Q92686
Solubility (25°C)	1600ng/ml
Sensitivity	2.42ng/ml
Entrez	4900
Expiration Date	Please enquire.