

Product Datasheet

Human Cathepsin G ELISA Kit (orb1213640)

Description	Human Cathepsin G, CTSG ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CTSG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,837
Uniprot ID	P08311
Solubility (25°C)	4800ng/L
Sensitivity	10.06ng/L
Entrez	1511
Expiration Date	Please enquire.