

Product Datasheet

Human 25-Hydroxyvitamin D-1 Alpha Hydroxylase ELISA Kit (orb1213644)

Description	Human 25-Hydroxyvitamin D-1 Alpha Hydroxylase, Mitochondrial, CYP27B1 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CYP27B1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,504
Uniprot ID	015528
Solubility (25°C)	4000ng/L
Sensitivity	9.47ng/L
Entrez	1594
Expiration Date	Please enquire.