

Product Datasheet

Human Regucalcin ELISA Kit (orb1213646)

Description	Human Regucalcin, RGN ELISA Kit
Reactivity	Human
Range	5-1000ng/ml
Tested Applications	ELISA
Target	RGN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,253
Uniprot ID	Q15493
Solubility (25°C)	1280ng/ml
Sensitivity	3.18ng/ml
Entrez	9104
Expiration Date	Please enquire.