

Product Datasheet

Human Rac-alpha Serine ELISA Kit (orb1213655)

Description	Human Rac-alpha Serine, Threonine-protein Kinase, AKT1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	AKT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,686
Uniprot ID	P31749
Solubility (25°C)	25.6ng/ml
Sensitivity	0.024ng/ml
Entrez	207
Expiration Date	Please enquire.