

## Product Datasheet

### Human Endothelial Protein C Receptor ELISA Kit (orb1213656)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Human Endothelial Protein C Receptor, EPCR ELISA Kit |
| <b>Reactivity</b>          | Human                                                |
| <b>Range</b>               | 0.5-100ng/ml                                         |
| <b>Tested Applications</b> | ELISA                                                |
| <b>Target</b>              | PROCR                                                |
| <b>Storage</b>             | 2-8°C                                                |
| <b>Note</b>                | For research use only                                |
| <b>Application notes</b>   | Sandwich                                             |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates               |
| <b>MW</b>                  | 26,671                                               |
| <b>Uniprot ID</b>          | <b>Q9UNN8</b>                                        |
| <b>Solubility (25°C)</b>   | 128ng/ml                                             |
| <b>Sensitivity</b>         | 0.27ng/ml                                            |
| <b>Entrez</b>              | <b>10544</b>                                         |
| <b>Expiration Date</b>     | Please enquire.                                      |