

Product Datasheet

Human Cerebellin-1 ELISA Kit (orb1213662)

Description	Human Cerebellin-1, CBLN1 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	CBLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,097
Uniprot ID	P23435
Solubility (25°C)	1600ng/L
Sensitivity	2.34ng/L
Entrez	869
Expiration Date	Please enquire.