

Product Datasheet

Human Protein S100-A10 ELISA Kit (orb1213664)

Description	Human Protein S100-A10, S100A10 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	S100A10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,203
Uniprot ID	P60903
Solubility (25°C)	32ng/ml
Sensitivity	0.027ng/ml
Entrez	6281
Expiration Date	Please enquire.