

Product Datasheet

Human Glycocalicin ELISA Kit (orb1213667)

Description	Human Glycocalicin, GC ELISA Kit
Reactivity	Human
Range	0.05-20ug/mL
Tested Applications	ELISA
Target	PAEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,624
Uniprot ID	P09466
Solubility (25°C)	24µg/mL
Sensitivity	0.028µg/mL
Entrez	5047
Expiration Date	Please enquire.