

Product Datasheet

Human 14-3-3 Protein Gamma ELISA Kit (orb1213676)

Description	Human 14-3-3 Protein Gamma, YWHAG ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	YWHAG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,303
Uniprot ID	P61981
Solubility (25°C)	96ng/ml
Sensitivity	0.16ng/ml
Entrez	7532
Expiration Date	Please enquire.