

Product Datasheet

Human Bcl-2-like Protein 11 ELISA Kit (orb1213684)

Description	Human Bcl-2-like Protein 11, BCL2L11 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	BCL2L11
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,171
Uniprot ID	043521
Solubility (25°C)	3200ng/L
Sensitivity	6.98ng/L
Entrez	10018
Expiration Date	Please enquire.