

Product Datasheet

Human Myelin Protein P0 ELISA Kit (orb1213685)

Description	Human Myelin Protein P0, MPZ ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	MPZ
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,555
Uniprot ID	P25189
Solubility (25°C)	4800ng/L
Sensitivity	9.14ng/L
Entrez	4359
Expiration Date	Please enquire.