

Product Datasheet

Human 14-3-3 Protein Epsilon ELISA Kit (orb1213686)

Description	Human 14-3-3 Protein Epsilon, YWHAE ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	YWHAE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,174
Uniprot ID	P62258
Solubility (25°C)	9600ng/L
Sensitivity	14.77ng/L
Entrez	7531
Expiration Date	Please enquire.