

Product Datasheet

Human Oxyntomodulin ELISA Kit (orb1213711)

Description	Human Oxyntomodulin, OXM ELISA Kit
Reactivity	Human
Range	5-1000ng/ml
Tested Applications	ELISA
Target	OMD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,492
Uniprot ID	Q99983
Solubility (25°C)	1280ng/ml
Sensitivity	2.41ng/ml
Entrez	4958
Expiration Date	Please enquire.