

Product Datasheet

Human Nischarin ELISA Kit (orb1213747)

Description	Human Nischarin, NISCH ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	NISCH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	166,629
Uniprot ID	Q9Y2I1
Solubility (25°C)	40ng/ml
Sensitivity	0.023ng/ml
Entrez	11188
Expiration Date	Please enquire.