

Product Datasheet

Human Acid Ceramidase ELISA Kit (orb1213750)

Description	Human Acid Ceramidase, ASAH1 ELISA Kit
Reactivity	Human
Range	0.5-180ng/ml
Tested Applications	ELISA
Target	ASAH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,660
Uniprot ID	Q13510
Solubility (25°C)	200ng/ml
Sensitivity	0.28ng/ml
Entrez	427
Expiration Date	Please enquire.