

Product Datasheet

Human Afamin ELISA Kit (orb1213751)

Description	Human Afamin, AFM ELISA Kit
Reactivity	Human
Range	5-600ng/ml
Tested Applications	ELISA
Target	AFAMIN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,069
Uniprot ID	P43652
Solubility (25°C)	640ng/ml
Sensitivity	2.38ng/ml
Entrez	173
Expiration Date	Please enquire.