

Product Datasheet

Human Sulphydryloxidase 1 ELISA Kit (orb1213754)

Description	Human Sulphydryloxidase 1, QSOX1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	QSOX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,578
Uniprot ID	000391
Solubility (25°C)	2400ng/L
Sensitivity	5.15ng/L
Entrez	5768
Expiration Date	Please enquire.