

Product Datasheet

Human Eosinophil Peroxidase ELISA Kit (orb1213757)

Description	Human Eosinophil Peroxidase, EPO ELISA Kit
Reactivity	Human
Range	0.2-60mU/ml
Tested Applications	ELISA
Target	EPX
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	81,040
Uniprot ID	P11678
Solubility (25°C)	80mU/ml
Sensitivity	0.12mU/ml
Entrez	8288
Expiration Date	Please enquire.