

Product Datasheet

Human Pro-Opiomelanocortin ELISA Kit (orb1213766)

Description	Human Pro-Opiomelanocortin, POMC ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	POMC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,424
Uniprot ID	P01189
Solubility (25°C)	9600ng/L
Sensitivity	14.84ng/L
Entrez	5443
Expiration Date	Please enquire.