

Product Datasheet

Human 14-3-3 Protein Eta ELISA Kit (orb1213772)

Description	Human 14-3-3 Protein Eta, YWHAH ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	YWHAH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,219
Uniprot ID	Q04917
Solubility (25°C)	96ng/ml
Sensitivity	0.14ng/ml
Entrez	7533
Expiration Date	Please enquire.