

Product Datasheet

Human Aldo-keto Reductase Family 1 ELISA Kit (orb1213774)

Description	Human Aldo-keto Reductase Family 1, Member B10, AKR1B10 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	AKR1B10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,020
Uniprot ID	O60218
Solubility (25°C)	12.8ng/ml
Sensitivity	0.021ng/ml
Entrez	57016
Expiration Date	Please enquire.