

Product Datasheet

Human Aldo-keto Reductase Family 1 ELISA Kit (orb1213775)

Description	Human Aldo-keto Reductase Family 1, Member B1, AKR1B1 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	AKR1B1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,853
Uniprot ID	P15121
Solubility (25°C)	40ng/ml
Sensitivity	0.02ng/ml
Entrez	231
Expiration Date	Please enquire.