

Product Datasheet

Human Kallistatin ELISA Kit (orb1213778)

Description	Human Kallistatin, SERPINA4 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	SERPINA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,542
Uniprot ID	P29622
Solubility (25°C)	24ng/ml
Sensitivity	0.022ng/ml
Entrez	5267
Expiration Date	Please enquire.