

Product Datasheet

Human Macrophage Scavenger Receptor 1 ELISA Kit (orb1213781)

Description	Human Macrophage Scavenger Receptor 1, MSR1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	MRC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	166,012
Uniprot ID	P22897
Solubility (25°C)	48ng/ml
Sensitivity	0.053ng/ml
Entrez	4360
Expiration Date	Please enquire.