

Product Datasheet

Human Alpha-methylacyl-CoA Recemase ELISA Kit (orb1213782)

Description	Human Alpha-methylacyl-CoA Recemase, AMACR ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	AMACR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,387
Uniprot ID	Q9UHK6
Solubility (25°C)	4800ng/L
Sensitivity	10.16ng/L
Entrez	23600
Expiration Date	Please enquire.