

Product Datasheet

Human Aromatase ELISA Kit (orb1213793)

Description	Human Aromatase, ARO ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	CYP19A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,883
Uniprot ID	P11511
Solubility (25°C)	64ng/ml
Sensitivity	0.11ng/ml
Entrez	1588
Expiration Date	Please enquire.